



HealthVision AI

AI-Based Visual Health Screening

Final Year Project Presentation

AI-Based Visual Health Screening System

A comprehensive AI-powered platform for non-invasive detection of Jaundice, Skin Diseases, and Nail Disorders using deep learning

Institution

Government Engineering College, Munger

Department of Computer Science & Engineering
(AI)

Project Supervisor

Dr. Saurabh Suman

Assistant Professor, CSE Department

Semester

7th Semester

B.Tech CSE (Artificial Intelligence)

Project Team Members

G

Gaurav Kumar
23151144901

N

Nitesh Kumar
22151144040

R

Rupesh Kumar
22151144010

I

Indrajeet Kumar
23151144906

Project Overview

A comprehensive AI-based health screening platform for non-invasive disease detection



Our Mission

Build a non-invasive, AI-powered health screening platform that uses deep learning to analyze images of eyes, face skin, and nails to detect potential health conditions — **without blood tests or lab visits**.



Eye Detection

Jaundice screening through sclera analysis



Face Detection

5-class skin disease classification



Nail Detection

17+ nail disorder identification

Key Statistics



Detection Modules

3



Detectable Conditions

24+



Analysis Time

<2s



Input Resolution

224×224



CNN Architectures

2



Powered by MobileNetV2 & EfficientNet Transfer Learning

Why This Project?

Addressing critical gaps in accessible health screening



Limited Rural Access

Access to **dermatological and ophthalmological screening** is severely limited in rural and underserved areas. Many patients must travel long distances just for basic health assessments, creating barriers to early detection.


⚠️ Rural healthcare gap affects millions



Invasive Tests Required

Jaundice detection traditionally requires blood tests to measure bilirubin levels. Our system enables **non-invasive screening** through simple sclera image analysis — no needles, no lab visits needed.


✅ Non-invasive alternative available



Specialist Shortage

Common skin conditions like **Acne, Eczema, Herpes, and Rosacea** require dermatological expertise. AI can assist in preliminary classification where specialists aren't available, bridging the expertise gap.

🧠 AI bridges expertise gap



Early Detection Gap

Nail abnormalities can indicate systemic health issues like liver disease, anemia, and autoimmune conditions. Early detection through image analysis prompts timely medical consultation, potentially saving lives.

💚 Early detection saves lives

Project Objectives

The key goals we set out to achieve with HealthVision AI

Core Objectives


- 1 Develop a multi-module AI platform capable of analyzing eye, face, and nail images for disease detection
- 2 Implement transfer learning using MobileNetV2 and EfficientNetB0 for high-accuracy classification
- 3 Build an intelligent eye detection pipeline using OpenCV Haar cascades with dual-strategy detection
- 4 Create a responsive web interface supporting both image upload and real-time camera capture
- 5 Ensure privacy by processing all data locally — no images stored or transmitted to third-party services
- 6 Provide detailed disease information alongside predictions for educational and screening purposes

Key Deliverables

-  **Three Trained Models**
.h5 format for eye, face, and nail detection
-  **Flask REST API Backend**
Prediction endpoints for each module
-  **Responsive Web Application**
Camera and upload support with SPA design
-  **Disease Information Database**
Comprehensive info on 20+ conditions
-  **Image Enhancement Pipeline**
Denoising and contrast normalization
-  **Complete Documentation**
Project report and deployment instructions

How It Works

A four-step pipeline from image capture to clinical prediction




Processing Speed
< 2 Seconds



Privacy
100% Local



Accuracy
Up to 90.62%



Device Support
All Devices

Detection Module 1: Eye / Jaundice Detection

Binary classification for jaundice screening through sclera analysis



About Jaundice

Jaundice is caused by elevated bilirubin levels, which turn the sclera (white of the eye) yellow. Our model analyzes eye images to detect this discoloration for early screening.

- ✓ Non-invasive detection method
- ✓ No blood tests required
- ✓ Early warning system



Technical Specifications



Architecture

MobileNetV2



Activation

Binary Sigmoid



Detection

Haar Cascade



Input Size

224×224



Model Parameters

~2.4 Million Parameters

Detection Pipeline



Upload Image



Eye Detection



Sclera Crop



AI Analysis



Result



Detection Module 2: Face / Skin Disease Detection

5-class classifier for facial skin conditions using EfficientNetB0

5 Skin Conditions Detected



Acne
Mild



Eczema
Moderate



Herpes
Moderate



Panu
Mild



Rosacea
Moderate

AI IN DERMATOLOGY

FASTER & MORE ACCURATE
SKIN DISEASE DETECTION



Technical Specifications

Architecture

EfficientNetB0

ImageNet pre-trained

Classification

5-Class Softmax

Multi-class output

Normalization

BatchNorm

Stabilized training

Dense Layers

256 + Dropout 0.5

Regularization

Model Performance

Validation Accuracy

72.15%

Training Images: 1,196

Validation Images: 298

Loss Function: Categorical Cross-Entropy

Detection Module 3: Nail Disease Detection

3-class nail disorder classifier for early health indicator detection

Why Nail Analysis?

Nail abnormalities can indicate systemic health issues including liver disease, anemia, fungal infections, and autoimmune conditions. Early detection through image analysis enables timely medical consultation.

- ✓ 17+ detectable conditions
- ✓ Systemic health indicators
- ✓ Non-invasive screening



3 Main Classes Detected



Healthy
Normal nail condition



Onychomycosis
Fungal infection



Psoriasis
Autoimmune condition

Technical Specifications

Architecture

MobileNetV2

Classification

3-Class Softmax

Input Size

224×224 pixels

Accuracy

85.71%

 Training: 933 images | Validation: 231 images

Technology Stack

Built with industry-standard tools and frameworks



Deep Learning



TensorFlow / Keras

Primary ML framework



MobileNetV2

Lightweight CNN architecture



EfficientNetB0

High-efficiency model



ImageNet Transfer Learning

Pre-trained weights



Web Backend



Python 3.10+

Programming language



Flask

Web framework



Flask-CORS

Cross-origin support



REST API / JSON

Data exchange format



Computer Vision



OpenCV

Image processing



Haar Cascades

Object detection



Pillow (PIL)

Image manipulation



NumPy

Numerical computing



Frontend



HTML5

Structure



CSS3

Styling



JavaScript ES6+

Interactivity



WebRTC

Camera access

Training Datasets

Curated medical image datasets for each detection module



Jaundice Eye

ezyZip Dataset

Classes 2

Jaundice, Normal

Training Images 136

Validation Images 32

Preprocessing

Rescale to [0, 1] ($\div$ 255)

 Split: 80% Train / 20% Val



Face Skin Disease

FaceZip Dataset

Classes 5

Acne, Eczema, Herpes, Panu, Rosacea

Training Images 1,196

Validation Images 298

Preprocessing

EfficientNet preprocess_input

 Split: 80% Train / 20% Val



Nail Disease

nailZip Dataset

Classes 3

Healthy, Onychomycosis, Psoriasis

Training Images 933

Validation Images 231

Preprocessing

Rescale to [0, 1] ($\div$ 255)

 Split: 80% Train / 20% Val



Data Augmentation (Training)



Rotation

20°



Shift

15%



Zoom

20%



Flip

Horizontal



Brightness

0.7-1.3



Shear

10%



Fill Mode

Nearest

Training Results & Accuracy

Two-stage training strategy with comprehensive performance metrics

**Jaundice Model**
MobileNetV2

Validation Accuracy
90.62%

 Loss Function
Binary Cross-Entropy

Train
136

Val
32

**Face Model**
EfficientNetB0

Validation Accuracy
72.15%

 Loss Function
Categorical Cross-Entropy

Train
1,196

Val
298

**Nail Model**
MobileNetV2

Validation Accuracy
85.71%

 Loss Function
Categorical Cross-Entropy

Train
933

Val
231

**Two-Stage Training Strategy**

1

Stage 1: Feature Extraction
Base frozen, LR: 1e-3, 15 epochs

2

Stage 2: Fine-Tuning
Top layers unfrozen, LR: 1e-5, 25 epochs

**Training Callbacks**


**EarlyStopping**
Patience: 5 epochs

**ReduceLROnPlateau**
Patience: 3, Factor: 0.3

**ModelCheckpoint**
Save best only

Application Flowchart

End-to-end data flow from user input to AI prediction output



Key Features & Capabilities

Six core capabilities that make HealthVision AI stand out



Smart Eye Detection

Haar cascade algorithm with dual-strategy detection: face-first and direct eye scanning with 35% padding for optimal sclera extraction.

- ✓ Intelligent fallback mechanism



Real-Time Camera

WebRTC camera integration with ROI guide overlay for precise image capture. Works seamlessly on both desktop and mobile browsers.

- ✓ Cross-platform compatibility



Deep Learning Models

MobileNetV2 and EfficientNetB0 transfer learning architectures pre-trained on ImageNet, fine-tuned for medical classification.

- ✓ State-of-the-art accuracy



Instant Results

Confidence scores, top-k predictions, and detected region previews generated in **under 2 seconds** per analysis.

- ✓ Lightning-fast inference



Privacy First

All processing happens locally on the server. No images are stored, shared, or transmitted to third-party services.

- ✓ 100% data security



Fully Responsive

Works seamlessly on desktops, tablets, and mobile phones with an adaptive UI that fits any screen size.

- ✓ Mobile-first design

Disease Information Database

Comprehensive coverage of 20+ detectable health conditions



Eye Conditions

2 Classes

Jaundice

Moderate

Elevated bilirubin causing yellow sclera

Normal

Healthy

Healthy eye with no discoloration



Skin Conditions

5 Classes

Acne

Clogged follicles

Mild

Eczema

Inflammatory patches

Moderate

Herpes

Viral blisters

Moderate

Panu

Fungal patches

Mild

Rosacea

Facial redness

Moderate



Nail Conditions

17+ Classes

Healthy

Onychomycosis

Psoriasis

Alopecia Areata

Beau's Lines

Bluish Nail

Clubbing

Darier's Disease

Half & Half

Koilonychia

Leukonychia

Lindsay's Nails

Muehrcke's

Onychogryphosis

Onycholysis

Pale Nail



Total Conditions

24+



Detailed Info

100%



Educational

Purpose



Thank You!

We appreciate your attention and interest in our project



Institution

**Government Engineering College,
Munger**

Department of CSE (AI)



Project Supervisor

Dr. Saurabh Suman
Assistant Professor



Project Team

Gaurav, Nitesh, Rupesh, Indrajeet
7th Semester CSE (AI)

Questions & Discussion

We would be happy to answer any questions about our project